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ATTENTION! Single Use Product. Do not use if the package is opened or damaged. Gamma radiation sterilized. Do not re-sterilize. For further safety, carefully read all directions herein before using this product.

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Implant Cylindrical UNIII HI features a threaded mechanically treated body. The threads (a) enhance implant placement and anchorage in the bone bed. The blasted implant body surface (b) demonstrates increased roughness but also increased surface area and anchorage, and increased bone integration. The collar shows (c) a finishing with low titanium roughness to improve its adherence to bone tissue. The three (d) uniformly distributed (120º) cutting edges at the apex improve implant placement. Available platform diameters (g): 3.50 mm for the 3.30 mm implant; 4.00 mm, for the 3.75 mm and 4.30 mm implant; 5.00 mm, for the

4.75 implant. The collar is barely divergent and the tapered apex is self-tapping (e). Internal hex makes the system easy to handle. The longitudinal axis of the upper part of the implant shows central threads (f) for the sitting of the prosthetic screw.

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Implant Cylindrical UNIII HE features a threaded body which (a) help enhance implant placement and anchorage in the bone bed. The blasted implant body surface (b) demonstrates increased roughness but also increased surface area and anchorage, increasing bone integration surface. The collar shows (c) a finishing with low titanium roughness to improve its adherence to bone tissue. The three (d) uniformly distributed (120°) cutting edges at the apex improve implant placement. Available platform diameters:

(g) 3.50 mm for the 3.30 mm implant; 4.00 mm, for the 3.75 mm and 4.30 mm implant; 5.00 mm, for the 4.75 mm implant. The collar is barely divergent and the tapered apex is self-tapping (e). Internal hex provides optimal prosthetic sitting and has been worldwide accepted for over 40 years. The longitudinal axis of the upper part of the implant shows central threads (f) for the sitting of the abutment screw. UN III HE Implant prosthetic sitting has been the preferred choice among professionals for decades worldwide.

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Implant Tapered UNIII HI features a threaded tapered body (a) which enhances implant placement and mechanical bone interlock. The implant body (b) is mechanically treated (blasting), which increases its surface and anchorage, and increases bone integration. The neck is barely divergent, and (c) shows a low titanium roughness to improve its adherence to bone tissue. The upper part of the implant shows micro-threads (d) which improve the load distribution on the cortical bone. The three cutting edges uniformly distributed (120°) at the apex (e) improve implant placement. Internal hex provides for optimal prosthetic sitting and has been worldwide accepted for over 40 years. The longitudinal axis of the upper part of the implant shows internal threads (f) for the sitting of the abutment screw. The tapered design of this implant favors the load distribution, and is especially indicated for low-density bone, or even for immediate loading.

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B. Medidas/Measures/Misure:

Diâmetros/Diameters/Diâmetros/Diametri: 3,5mm; 4,0mm; 5,0mm.

Comprimentos/Length/Largo/Lunghezze: 5mm; 6mm; 7mm; 8mm; 9mm; 10mm; 11mm; 12mm; 13mm; 14mm; 15mm; 16mm; 17mm; 18mm

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Implant Tapered UNIII HE – External Hexagon - features a threaded tapered body (a) which enhances the implant placement and bone interlock. The implant body (b) is mechanically treated (blasting), which increases its surface and enhances bone integration. The neck is barely divergent, and (c) shows a low titanium roughness to improve its adherence to bone tissue. The upper part of the implant body shows micro-threads (d) which improve the load distribution on the cortical bone. The three cutting edges uniformly distributed (120°) at the apex (e) improve implant placement. The tapered design of this implant favors the load distribution and is especially indicated for low-density bone or even for immediate loading. The HE UN III Tapered The HE UN III Tapered Implant combines both the efficiency of the external hex, which has worldwide recognition, and the user friendly tapered design.

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A) Product description The Implant Morse Taper UNIII AR 'Due Cone' features a threaded tapered body (a) which enhances implant placement and mechanical bone interlock. The implant body is mechanically treated (blasting), which increases its surface and anchorage, and increases bone integration. The surface treatment includes the entire neck, since this implant has to be placed 1 to 2 mm below the bone crest. The three uniformly distributed cutting edges at the apex (120°) improve implant placement. The neck (c) is barely divergent. The Morse Taper connection favors its use for single units. It has an internal hexagon (e) of 2.00 mm used for placement and as indexer for prosthetic units. The longitudinal axis of the upper part of the implant shows central

threads (d) for sitting of the abutment screw. The design of this implant favors the load distribution and is especially indicated for low-density bone. It is machined from commercially pure titanium according to ASTM F-67.

B) Accompanying Parts and Accessories

The product is packed individually Each pack contains: - one Implant Morse Taper UNIII AR 'Due Cone' - one cover screw

C) Technical specifications and features

Placement Torque: 60Ncm

Dimensions Ø of the base (mm): 3,5; 4,0; 4,5; 5,0; 5,5

Length (mm): 5; 6; 7; 8; 9; 10; 11; 12; 13; 14; 15; 16; 17; 18

Raw Material: The raw material used for the "Implant Morse Taper UNIII AR 'Due Cone'" is Titanium according to ASTM F-67.

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Special product transportation, storage and maintenance conditions:
Storage: Keep in site protected from humidity and sun, and in its original package.
Transportation: When transporting, protect package from humidity.
Maintenance: Keep in site protected from humidity and sun.

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Warnings and / or precautions:

For the safe and effective use of Implacil[®] implants, we strongly recommend specific training. Highly specialized and complex surgical techniques are required to place such dental implants. An improper technique may cause implant failure and / or loss of supporting bone or any other kind of negative effects. Implant patients should be thoroughly selected. Age and the presence of parafunctional habits should be strongly analyzed. Visual inspection and periapical X-rays are critical to check anatomy, occlusion, periodontal and bone status. The overall patient health should be adequate for the placement of implants. Do not sterilize or re-use Implacil[®] implants, the risk of reuse is biological contamination! Do not use an open or damaged package. The temperature of 47°C applied for one minute might cause bone necrosis. At the temperature of 50°C for more than one minute the bone around the implant reaches the point of non-regeneration. The use of progressive diameter cutting burs, adequate speed (1500 to 2000 RPM/min), and abundant irrigation prevents excessive heating. This product is contraindicated for patients exhibiting signs of allergy or hypersensitivity to the chemical ingredients of the material: titanium.

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

Product usage directions:

- Open the outside package. Check whether the blister is perfectly sealed. Do not use if the package is damaged.
- Start the sequence of drills to prepare the implant site. Open the blister right before implant insertion. C) Remove the plastic capsule from the blister and open it carefully. Do not discard the capsule lid. The cover screw is inside it.
- Remove the implant from the capsule using the driver tool.
- Insert the implant in the recently drilled site. Maximum insertion speed is 20 RPM. If manually done, maximum torque is 60N.mm.
- After inserted, remove the mount using 1.17 mm hexagonal tool.
- Mount being removed.
- Remove the cover screw from the capsule lid.
- Seal with the cover screw.
- Installation complete.

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Medical device safety and efficacy requirements, and possible undesirable secondary effects.

5.1 – Indication, function or usage of the product. Indication: Patients needing total or partial oral rehabilitation.

Function: replace the root of the lost tooth.

5.2 – Undesirable secondary or side effects and contraindications: Implant surgery is a highly specialized and complex procedure and should only be performed by professionals with appropriate technical qualification. Overloading may cause excessive bone loss or implant fracture. Physiology and anatomy may impair dental implants performance. These aspects should be carefully considered in patients with the following conditions: low quality bone; poor oral hygiene; medical conditions such as blood diseases or uncontrollable hormone disorders; metabolic or systemic diseases affecting bone or wound healing, uncontrolled diabetes and coagulation disorders; chemo or radiation therapy; untreatable occlusion or articulation problems; insufficient bone height or width; tobacco or alcohol abuse. Potential negative events associated to dental implants may include: non integration; loss of integration; dehiscence needing bone graft; perforation of maxillary sinus, lingual face or lower alveolar canal; persistent pain, numbness, paresthesia; excessive bone loss needing intervention; and implant fracture.

5.3 – Product safety and efficacy: Universal III implant is totally safe, once safety rules are followed in all procedures. The product was developed with medical, non-toxic material, made from ASTM F-67 titanium. When stored according to item 2, the product will not loose or change its characteristics and dimensions. The product is sterile and adequate for use for up to 2 years as from the sterilization date.

Magnetic Resonance Imaging (MRI): implacil De Bortoli did not evaluate the safety of its dental implants with the magnetic resonance environment. The safety in relation to the magnetically induced displacement force and the extent to which radiofrequency (RF) induced heating is able to affect patient safety have not been evaluated. The safety of Implacil De Bortoli dental implants in the MRI environment is unknown. Performing an MRI on a patient with this device can result in damage to the patient.

Implacil De Bortoli recommends that medical professionals and the manufacturer of the MRI equipment be consulted for compatibility with implants before use. We inform that Implacil implants are machined in Titanium ASTM F67 Gr.4.

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Installation or Connection to Other Products: To be implanted, the UNIVERSAL III IMPLANT has to be connected to an adequate wrench. After integration the prosthetic components are connected to obtain patient's oral restoration.

Attention: The use of any part, accessory or material not specified or provided here is of total responsibility of the user.

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Discarding Precautions: To prevent environmental contamination or the inadequate use of the Prosthetic Abutments, when it is rendered useless, it should be adequately discarded in a site for potentially contaminated materials.

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
Warranty Term: IMPLACIL DE BORTOLI - MATERIAL ODONTOLOGICO S.A. assures to the customer of this product, warranty against any material or manufacturing defect for a period of one year. The warranty will be lost in the case of the original customer resale the product. This product has undetermined validity when adequately transported, stored and maintained according to these Usage Directions.

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Instructions for Use for all countries except Brazil

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