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IFU-017 INSTRUCTIONS FOR USE/ IFU-017 INSTRUCCIONES DE USO/ IFU-017 ISTRUZIONI PER L'USO

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<p>4350 - Transfer Abutment Closed Tray HE 5,0 204675 – Transfer HE Cubeta Cerrada 3,5 4336 – Transfer HE Cubeta Cerrada 4,0 4350 – Transfer HE Cubeta Cerrada 5,0 204675 – Transfer HE Moldeador Cerrado 3,5 4336 – Transfer HE Moldeador Cerrado 4,0 4350 - Transfer HE Moldeador Cerrado 5,0</p>		<p>14854 - Transfer HI Porta Impronta Aperto 3,5 15899 - Transfer HI Porta Impronta Aperto 4,0 14779 - Transfer HI Porta Impronta Aperto 5,0</p>	
<p>4404 - Transfer HI 3,5 4411 - Transfer HI 4,0 4428 - Transfer HI 5,0/ 4404 - Transfer Abutment Closed Tray HI 3,5 4411 - Transfer Abutment Closed Tray HI 4.0 4428 - Transfer Abutment Closed Tray HI 5.0 4404 - Transfer HI Cubeta Cerrada 3,5 4411 - Transfer HI Cubeta Cerrada 4,0 4428 - Transfer HI Cubeta Cerrada 5,0 4404 - Transfer HI Moldeador Cerrado 3,5 4411 - Transfer HI Moldeador Cerrado 4,0 4428 - Transfer HI Moldeador Cerrado 5,0</p>		<p>13512 - Transfer Mini Cônico Moldeira Aberta/13512 - Mini Tapered Open Tray Transfer/13512 - Transfer Mini Cônico Cubeta Abierta/13512 - Transfer Mini Conico Porta Impronta Aperto</p>	
<p>4282 - Transfer Mini Cônico Moldeira Fechada/ 4282 - Mini Tapered Transfer Closed Tray /4282 - Transfer Mini Cônico Cubeta Cerrada /4282 - Transfer Mini Conico Porta Impronta Chiuso</p>		<p>4176 - Transfer Pilar Hexagonal/4176 - Hexagonal Abutment Transfer/4176 - Transfer Pilar Hexagonal/4176 - Transfer Pilastro Esagonale</p>	
<p>214223 - Transfer HE Moldeira Aberta 4.0 ST/214223 - HE Open Tray Transfer 4.0 ST/214223 - Transfer HE Cubeta Abierta 4.0 ST/214223 - Transfer HE Porta Impronta Aperto 4.0 ST</p>		<p>214254 - Transfer HE Moldeira Fechada 4.0 ST/214254 - HE Closed Tray Transfer 4.0 ST/214254 - Transfer HE Cubeta Cerrada 4.0 ST/214254 - Transfer HE Porta Impronta Chiuso 4.0 ST</p>	
<p>214247 - Transfer HI Moldeira Aberta 4.0 ST/214247 - HI Open Tray Transfer 4.0 ST/214247 - Transfer HI Porta Impronta Aperto 4.0 ST</p>		<p>214278 - Transfer HI Moldeira Fechada 4.0 ST/214278 - HI Closed Tray Transfer 4.0 ST/214278 - Transfer HI Porta Impronta Chiuso 4.0 ST/ HI Cubeta Cerrada 4.0 ST/</p>	

<p>228923 - Transfer Implante Cone Morse AR Moldeira Fechada/228923 - Morse Taper Abutment AR Closed Tray Transfer/228923 - Transfer Implante Cono Morse AR Cubeta Cerrada/228923 - Transfer Impianto Cono Morse AR Porta Impronta Chiuso</p>		<p>228930 - Transfer Implante Cone Morse AR Moldeira Aberta/228930 - Morse Taper Abutment AR Open Tray Transfer/228930 - Transfer Implante Cono Morse AR Cubeta Abierta/228930 - Transfer Impianto Cono Morse AR Porta Impronta Aperto</p>	
<p>21711 - Transfer Clássico 3.25 21713 - Transfer Clássico 3.75 217156 - Transfer Clássico 4.25/ 21711 - Classic Transfer 3.25 21713 - Classic Transfer 3.75 217156 - Classic Transfer 4.25/ 21711 - Transfer Clásico 3.25 21713 - Transfer Clásico 3.75 217156 - Transfer Clásico 4.25/ 21711 - Transfer Classico 3.25 21713 - Transfer Classico 3.75 217156 - Transfer Classico 4.25</p>		<p>25466 - Transfer Micro Cônico Moldeira Aberta 25466 - Micro Tapered Open Tray Transfer 25466 - Transfer Micro Cônico Cubeta Abierta 25466 - Transfer Micro Conico Porta Impronta Aperto</p>	
<p>26464 - Transfer Micro Cônico Moldeira Fechada/26464 - Micro Tapered Closed Tray Transfer /26464 - Transfer Micro Cônico Cubeta Cerrada/26464 - Transfer Micro Conico Porta Impronta Chiuso</p>		<p>231497 - Transfer Implante Cone Morse ST 231497 - Transfer Implant Cone Morse ST 231497 - Transfer Implante Cono Morse ST 231497 - Transfer Impianto Cono Morse ST</p>	

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EN 1- COMMERCIAL PRESENTATION FORMS OF THE PRODUCT

The Impression Transfers are supplied in individual packs.

The primary packaging consists of a blister produced with virgin PET plastic, and it is sealed with Tyvek.

List of items contained inside the package

- 1 Transfer

- 3 Record Tags

- This instruction for use will be currently available in electronic format on the website of Implacil De Bortoli (www.debortoli.com.br/manuais.asp) according to the effective legislation on the availability of instruction for use in no-printed format of healthcare products – IN Normative Instruction No. 4 dated June 15, 2012.

- The printed instruction for use will be available when requested by a customer and sent free of charge when requested through Customer Service: 0800-779-0011.

To locate it, check the registration number of the product at ANVISA and the version of the instruction indicated on the Label of the product purchased.

- Adobe Acrobat Reader is necessary in order to read the instructions for use on the website.

The individual packs are square in shape with dimensions of 80mmx80mm.

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Transfer HI Moldeira Fechada: 3,5 / 4,0 / 5,0

Transfer Clássico: 3,25 / 3,75 / 4,25 mm

EN Commercial Models with dimensional variations:

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EN Other models do not have dimensional variations.

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EN 2- INDICATION OF USE AND PURPOSE

The Impression Transfer is indicated for oral rehabilitation cases where it is necessary to reproduce the positioning of the implant or prosthetic abutment in the patient's mouth that will be worked on in the clinic or laboratory. Each transfer is manufactured for a specific type of implant or prosthetic abutment.

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EN 3- MODE OF PRODUCT USE

- Previously select the Transfer suitable for the impression technique and the type of impression and component to be used.
- The Implacil Impression Transfer is packed inside a blister sealed with Tyvek.
- Check if the pack is in perfect conditions. Reject the product if the pack is damaged.
- Check if it is perfectly sealed.
- Remove the transfer from inside the blister.
- The Transfer is previously sanitized in a classified room. Follow the stages below if you wish to sanitize again:

CLEANING

The cleaning of a surgical instrument that has been in contact with tissue must be conducted in such a way that all parts of the item are exposed to the cleaning process.

For cleaning Implacil surgical instruments use only the following manual cleaning process.

Automated cleaning has not been validated. Do not use automated cleaning for Implacil surgical instruments.

1. Initial rinse with cold water. The initial cold rinse is conducted in order to remove any gross contamination present on the item. Cold water is used as temperatures above 60°C/140°F will coagulate proteins, making it difficult to remove them from contaminated items.
2. Prepare the enzymatic detergent according to the recommendation*.
3. Soak the products in the prepared detergent solution and soak for at least 5 minutes. While soaking, brush items with a soft-bristle brush.
4. Rinse: after removing the instruments from the detergent solution, rinse thoroughly for at least 3 times (for at least 1 minute) in running water.
5. Visually inspect each product for visible dirt. If there is any visible dirt, repeat the cleaning process (step 3 and 4) until no more visible dirt is left.

*Recommendation for the cleaning step

Always use distilled, deionized or demineralized water for this procedure.

Never use saline solutions, especially sodium hypochlorite and saline solution, disinfectants, hydrogen peroxide or alcohol to clean or rinse the instruments.

Never use steel wool or sponges and abrasive products, so that the instruments are not damaged.

Cleaning: Enzymatic Detergent Composition: Protease, Lipase and Amylase enzymes

Dilution: 1mL for each liter of water;

Action time: 5 min (Manual);

Water temperature: 40°C to 55°C (Manual);

Note: The cleaning process must be conducted in such a way that all parts of the surgical instrument are exposed. Contamination remaining on an item after cleaning reduces the effectiveness of any subsequent decontamination or sterilization process.

Clean instruments as soon as possible after use. Do not allow blood and debris to dry on the instruments. If cleaning needs to be delayed, place groups of instruments in a covered container with an enzyme solution or appropriate detergent to delay drying.

Surgical instruments must be replaced as soon as the professional perceives their reduced effectiveness.

STERILIZATION

• The product is not sterile, and this procedure should be carried out by the professional according to the following stages:

This product must be properly sanitized and sterilized before use. Please use only steam sterilization according to the following parameters for sterilization:

Sterilization Parameters	
Temperature	132°C
Sterilization time	20 min
Drying time	20 - 60 min

Do NOT open the autoclave to avoid rapid condensation, i.e., allow all the steam to come out first in order to allow the Drying Cycle to be completed on its own. Clean the autoclave thoroughly and periodically.

- Thread the Impression Transfer on the component that you wish to transfer to the cast.
- Use the impression technique to obtain an impression with the exact position of the Impression Transfer.
- Remove the Impression Transfer from the patient's mouth and insert it in the position in the impression obtained.
- Insert the analog of the component to be transferred in the impression threading it in the Impression Transfer.
- Pour out the plaster to obtain the working cast containing the replica of the component in the position that this is originally found in the mouth.

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EN 4- COMPOSITION: Stainless Steel 303 (ASTM F899) is used as raw material for the Transfers

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EN 6- STORAGE CONDITIONS

Keep it in the original packaging protected from direct sunlight and rain.

The product in question is not thermosensitive and does not present any variation in environments with high or low humidity levels; therefore, it can be sterilized in an autoclave. Even so, it is recommended to store the product between 15° and 35° C and humidity between 45% and 80% in order to ensure the stability of the packaging.

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EN 7- TRANSPORT CONDITIONS

The Transfers should be transported in a clean and dry place and at room temperature. The transport should be performed appropriately to avoid falls and damages in their original pack.

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EN 8- MANIPULATION CONDITIONS

Handling should be performed with the proper precautions to avoid falls and damages, which must be performed on shelves and counters according to the activity that is being executed.

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EN 9- WARNINGS

The Implacil Impression Transfers should only be used by professionals specialized in the field. It is necessary to elaborate a prosthetic execution plan, observing the diameter, angle and length for adaptation of the prosthetic component to the implant in order to avoid future complications. The product is not sterile, thus this procedure must be performed by the professional before any intervention on the patient.

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EN 10- PRECAUTIONS

The Implacil Impression Transfers were designed for single use, not recommending reprocessing. Do not use it at all if the pack is damaged. The product is not sterile, thus this procedure must be performed by the professional before any intervention on the patient.

The professional must inform the patient of the need for periodic monitoring, to avoid physical and mechanical stresses in order to not subject the product to undue stresses, as well as the need for correct hygiene.

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EN 11 – CONTRAINDICATIONS: Not applicable.

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
EN 14 - Warranty Term

IMPLACIL DE BORTOLI - MATERIAL ODONTOLOGICO S.A. guarantees this product against any material or manufacturing defect that may appear within 365 days, counting from the original date of purchase.

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Mário Sérgio De Bortoli

Responsável Legal / Legal Representative / Responsable Legal / Responsabile Legale

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EN Authorized Representative in the European Community:

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