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As embalagens individuais possuem o formato retangular de 50mmx80mm.

EN 1- COMMERCIAL PRESENTATION FORMS OF THE PRODUCT

The Impression Transfers are supplied in individual packs.

The primary packaging consists of a blister produced with virgin PET plastic, and it is sealed with Tyvek.

List of items contained inside the package

- 1 Transfer

- 3 Traceability Tags

- This instruction for use will be currently available in electronic format on the website of Implacil De Bortoli (www.debortoli.com.br/manuais.asp) according to the effective legislation on the availability of instruction for use in no-printed format of healthcare products – Normative Instruction IN No. 4 dated June 15, 2012.

- The printed instruction for use will be available when requested by a customer and sent free of charge when requested through Customer Service: 0800-779-0011.

- To locate it, check the number of the instructions for use: "IFU-XXX" and the version of the instruction that is indicated on the Label of the purchased product.

- Adobe Acrobat Reader is necessary in order to read the instructions for use on the website.

The dimensional features of the plastic transfers of this process vary in diameter, X height, where the first measurement described corresponds to the diameter, which is 3.3 or 4.5 and the second measurement corresponds to the height that may be 4 or 6 mm.

The transfer models of this process are orange in color.

Dimensional features of the packaging:

The individual packs are rectangular in shape with dimensions of 50mmx80mm.

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EN 3- MODE OF PRODUCT USE

• Perform the prior selection of the Transfer suitable for the impression technique and the type of impression and component to be used.

- The Implacil Impression Transfer is packed inside a blister sealed with Tyvek.
- Check if the pack is in perfect conditions. Reject the product if the pack is damaged.
- Check if it is perfectly sealed.

• Remove the transfer from inside the blister.

• The Transfer is previously sanitized in a classified room. Follow the stages below if you wish to sanitize again:

- Initial rinse of the item with cold water.

- Wash with neutral to slightly basic detergent

- Final rinse with deionized or purified water by reverse osmosis.

• The product is not sterile, and this procedure should be carried out by the professional according to the following stages:

- Sterilization temperature: 132°C, Sterilization cycle: 20 minutes and Drying cycle: 20 – 60 minutes. The items should be arranged inside the container (perforated on the cover and lower part), or tray or in the box, or surgical envelop, provided they are arranged in a way that they can be prepared for the sterilization process. Do NOT open the autoclave to avoid rapid condensation, i.e., allow all the steam to come out first in order to allow the Drying Cycle to be completed on its own. Clean the autoclave thoroughly and periodically.

• Thread the Impression Transfer on the component that you wish to transfer to the cast.

• Use the impression technique to obtain an impression with the exact position of the Impression Transfer.

• Remove the Impression Transfer from the patient's mouth and insert it in the position in the impression obtained.

• Insert the analog of the component to be transferred in the impression threading it in the Impression Transfer.

• Pour out the plaster to obtain the working cast containing the replica of the component in the position that this is originally found in the mouth.

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EN 5- EXPIRY DATE: 5 years as from the date of manufacture.

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EN 6- STORAGE CONDITIONS

Keep it in the original packaging protected from direct sunlight and rain.

The product in question is not thermosensitive and does not present any variation in environments with high or low humidity levels; therefore, it can be sterilized in an autoclave. Even so, it is recommended to store the product between 15° and 35° C and humidity between 45% and 80% in order to ensure the stability of the packaging.

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EN 7- TRANSPORT CONDITIONS: The transfers should be transported in a clean and dry place and at room temperature. The transport should be performed appropriately to avoid falls and damages in their original pack.

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EN 8- MANIPULATION CONDITIONS: Handling should be performed with the proper precautions to avoid falls and damages, which must be performed on shelves and counters according to the activity that is being executed.

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EN 9- WARNINGS

The Implacil Impression Transfers should only be used by professionals specialized in the field. It is necessary to elaborate a prosthetic execution plan, observing the diameter, angle and length for adaptation of the prosthetic component to the implant in order to avoid future complications. The product is not sterile, thus this procedure must be performed by the professional before any intervention on the patient.

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EN 10- PRECAUTIONS

The Implacil Impression Transfers were designed for single use, not recommending reprocessing. Do not use it at all if the pack is damaged. The product is not sterile, thus this procedure must be performed by the professional before any intervention on the patient.

The professional must inform the patient of the need for periodic monitoring, to avoid physical and mechanical stresses in order to not subject the product to undue stresses, as well as the need for correct hygiene.

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EN 11- CONTRAINDICATIONS: Not applicable.

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
EN 14 – Warranty Term

IMPLACIL DE BORTOLI - MATERIAL ODONTOLOGICO S.A. guarantees this product against any material or manufacturing defect that may appear within 365 days, counting from the original date of purchase.

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EN Authorized Representative in the European Community:

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