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IFU-040 INSTRUCTIONS FOR USE / IFU-040 INSTRUCCIONES DE USO / IFU-040 ISTRUZIONI D'USO

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<p>30771 TRANSFER DIGITAL HE 3.5 AR EN TRANSFER DIGITAL HE 3.5 AR ES TRANSFER DIGITAL HE 3.5 AR ITA TRANSFER DIGITALE HE 3.5 AR</p>	 <p>HE35</p>
<p>30773 TRANSFER DIGITAL HE 4.0 AR EN TRANSFER DIGITAL HE 4.0 AR ES TRANSFER DIGITAL HE 4.0 AR ITA TRANSFER DIGITALE HE 4.0 AR</p>	 <p>HE40</p>
<p>30775 TRANSFER DIGITAL HE 5.0 AR EN TRANSFER DIGITAL HE 5.0 AR ES TRANSFER DIGITAL HE 5.0 AR ITA TRANSFER DIGITALE HE 5.0 AR</p>	 <p>HE50</p>
<p>30777 TRANSFER DIGITAL HI 3.5 AR EN TRANSFER DIGITAL HI 3.5 AR ES TRANSFER DIGITAL HI 3.5 AR ITA TRANSFER DIGITALE HI 3.5 AR</p>	 <p>HI35</p>
<p>30779 TRANSFER DIGITAL HI 4.0 AR EN TRANSFER DIGITAL HI 4.0 AR ES TRANSFER DIGITAL HI 4.0 AR ITA TRANSFER DIGITALE HI 4.0 AR</p>	 <p>HI40</p>
<p>30781 TRANSFER DIGITAL HI 5.0 AR EN TRANSFER DIGITAL HI 5.0 AR ES TRANSFER DIGITAL HI 5.0 AR ITA TRANSFER DIGITALE HI 5.0 AR</p>	 <p>HI50</p>
<p>30789 TRANSFER DIGITAL IDEALE 3.3X4 AR EN TRANSFER DIGITAL IDEALE 3.3X4 AR ES TRANSFER DIGITAL IDEALE 3.3X4 AR ITA TRANSFER DIGITALE IDEALE 3.3X4 AR</p>	 <p>ID334</p>

<p>30791 TRANSFER DIGITAL IDEALE 3.3X6 AR EN TRANSFER DIGITAL IDEALE 3.3X6 AR ES TRANSFER DIGITAL IDEALE 3.3X6 AR ITA TRANSFER DIGITALE IDEALE 3.3X6 AR</p>	 <p>ID336</p>
<p>30793 TRANSFER DIGITAL IDEALE 4.5X4 AR EN TRANSFER DIGITAL IDEALE 4.5X4 AR ES TRANSFER DIGITAL IDEALE 4.5X4 AR ITA TRANSFER DIGITALE IDEALE 4.5X4 AR</p>	 <p>ID454</p>
<p>30795 TRANSFER DIGITAL IDEALE 4.5X6 AR EN TRANSFER DIGITAL IDEALE 4.5X6 AR ES TRANSFER DIGITAL IDEALE 4.5X6 AR ITA TRANSFER DIGITALE IDEALE 4.5X6 AR</p>	 <p>ID456</p>
<p>30826 TRANSFER DIGITAL MICRO CONICO R EN TRANSFER DIGITAL MICRO CONICO R ES TRANSFER DIGITAL MICROCONÍCO R ITA TRANSFER DIGITALE MICRO CONICO R</p>	 <p>MC C</p>
<p>30783 TRANSFER DIGITAL MINI CONICO R EN TRANSFER DIGITAL MINI CONICO R ES TRANSFER DIGITAL MINICÓNICO R ITA TRANSFER DIGITALE MINI CONICO R</p>	 <p>MN C</p>
<p>31737 TRANSFER PILAR DIGITAL 3.8 AR EN TRANSFER PILAR DIGITAL 3.8 AR ES TRANSFER PILAR DIGITAL 3.8 AR ITA TRANSFER PILAR DIGITALE 3.8 AR</p>	 <p>PD38</p>
<p>31739 TRANSFER PILAR DIGITAL 4.5 AR EN TRANSFER PILAR DIGITAL 4.5 AR ES TRANSFER PILAR DIGITAL 4.5 AR ITA TRANSFER PILAR DIGITALE 4.5 AR</p>	 <p>PD45</p>
<p>30797 TRANSFER DIGITAL SMART 3.5X4 AR EN TRANSFER DIGITAL SMART 3.5X4 AR ES TRANSFER DIGITAL SMART 3.5X4 AR ITA TRANSFER DIGITALE SMART 3.5X4 AR</p>	

<p>30799 TRANSFER DIGITAL SMART 3.5X6 AR EN TRANSFER DIGITAL SMART 3.5X6 AR ES TRANSFER DIGITAL SMART 3.5X6 AR ITA TRANSFER DIGITAL SMART 3.5X6 AR</p>	
<p>30801 TRANSFER DIGITAL SMART 4.5X4 AR EN TRANSFER DIGITAL SMART 4.5X4 AR ES TRANSFER DIGITAL SMART 4.5X4 AR ITA TRANSFER DIGITALE SMART 4.5X4 AR</p>	
<p>30803 TRANSFER DIGITAL SMART 4.5X6 AR EN TRANSFER DIGITAL SMART 4.5X6 AR ES TRANSFER DIGITAL SMART 4.5X6 AR ITA TRANSFER DIGITALE SMART 4.5X6 AR</p>	
<p>34397 - TRANSFER DIGITAL PILAR STRETTO 3.0x4.8 EN TRANSFER DIGITAL PILAR STRETTO 3.0x4.8 ES TRANSFER DIGITAL PILAR STRETTO 3.0x4.8 ITA TRANSFER DIGITALE PILAR STRETTO 3.0x4.8</p>	
<p>34563 - TRANSFER DIGITAL PILAR STRETTO 3.0x3.5 EN TRANSFER DIGITAL PILAR STRETTO 3.0x3.5 ES TRANSFER DIGITAL PILAR STRETTO 3.0x3.5 ITA TRANSFER DIGITALE PILAR STRETTO 3.0x3.5</p>	

Registro ANVISA nº: 80259860043
EN ANVISA Registration No.: 80259860043
ES Registro ANVISA n.º: 80259860043
ITA Registrazione ANVISA n. 80259860043

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As embalagens individuais possuem o formato retangular de 50mmx80mm.

EN 1- FORMS OF PRODUCT BUSINESS PRESENTATION

The Digital Transfers are presented in individual packages.

The primary packaging consists of a blister pack produced from pure PET-type plastic, and is sealed with Tyvek.

List of items contained inside the package

- 01 Transfer
- 01 Screw (where applicable)
- 03 Traceability Tags

- This instruction for use will be available in electronic format on the website of Implacil De Bortoli (www.debortoli.com.br/manuais.asp) in accordance with the legislation in force on the availability of instructions for use in a non-printed format for health products - Normative Instruction IN No. 4 of June 15, 2012.

- The printed instruction of use will be available to the customer who requests it, and will be sent free of charge, when requested through the Customer Service: 0800-779-0011.

To find it, check the instruction number: "IFU-XXX" and the version of the instruction that is indicated on the Label of the purchased product.

- Adobe Acrobat Reader is required to read the instructions on the site.

The dimensional characteristics of the transfers of this process, vary in some models in diameter X height, where the first measure described corresponds to the diameter, and the second measure corresponds to the height. For the models that only have one measure, this one refers to the diameter. The unit of measurement is always mm.

Dimensional characteristics of packaging:

Individual packaging have the rectangular format of 50mmx80mm.

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EN 2- INDICATION OF USE AND PURPOSE

The Digital Transfer must be used during the model scanning in order to enable the alignment with the original library, indicating the correct position of the analog/implant.

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EN 3- PRODUCT USE MODE

- Pre-select the appropriate Transfers for the component to be used.
- The Implacil Digital Transfer is packaged inside a Tyvek-sealed blister.
- Check if the packaging is in perfect condition. Discard the product if its packaging is damaged.
- Check if it is perfectly sealed.
- Remove the Transfers from inside the blister.
- The Transfer is previously sanitized in a classified room. If you wish to sanitize again, here are the steps:
 - Initial rinse the item with cold water.
 - Washing with a neutral to slightly basic detergent
 - Final rinse with deionized or reverse osmosis purified water.
- The product is not sterile and this procedure must be performed by the professional following the steps:
 - Sterilization temperature: 132°C, Sterilization cycle: 20 minutes and Drying cycle: 20 – 60 minutes. The items should be arranged inside the container (perforated on the cover and on the bottom), or tray or in the box, or in the surgical envelope; as long as they are accommodated in such a way that they can be prepared for the sterilization process. DO NOT open the autoclave to prevent rapid condensation; i.e. let all steam escape first making the Drying Cycle complete by itself. Clean the autoclave rigorously and periodically.

Step 1: Definition of the CAD-CAM transfer to be used according to the characteristics of the implant or prosthetic intermediate;

Step 2: Screwing of the transfer over the implant or prosthetic intermediate; (where applicable).

Step 3: Scanning of the intraoral transfer or on the model.

Check the software for library compatibility with our .STL files or corresponding files in other planning software according to the codes described for CadCam use available in the instruction for use.

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EN 5 - EXPIRATION DATE: 5 years from the date of manufacture.

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EN 6 - STORAGE CONDITIONS

Keep in its original packaging in a place protected from direct sunlight and rain.

The product in question is not thermosensitive and does not present any variation in environments with high or low levels of humidity. Even so, it is recommended that the product be kept between 15° and 35° C and at a humidity between 45% and 80%, to guarantee the stability of the package.

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EN 7 - TRANSPORT CONDITIONS

Transfers must be transported in a clean, dry place at room temperature. Transportation must be done properly, to avoid dropping and damage to its original packaging.

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EN 8 - HANDLING CONDITIONS

Handling should be carried out with due care to avoid falls and damage, and it should be done on the shelves and counters, according to the activity being performed.

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EN 9 - WARNINGS

The use of Implacil Digital Transfers must be made only by specialized professionals in the area. It is necessary to elaborate a prosthetic execution plan, observing the diameter, angle and length for the adaptation of the prosthetic component to the implant or prosthetic intermediate in order to avoid future complications.

The product is not sterile, and this procedure should be performed by the professional before any intervention on the patient.

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EN 10 - PRECAUTIONS

Implacil's Digital Transfers were designed for single use, not being recommended for reprocessing.

Do not use it under any circumstances if the packaging is violated. The product is not sterile, and this procedure should be performed by the professional before any intervention on the patient.

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EN 11 – CONTRAINDICATIONS: Not applicable.

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EN 14 – Warranty Term

IMPLACIL DE BORTOLI - MATERIAL ODONTOLOGICO S.A. assures the owner of this product, warranty against any material or manufacturing defects that may appear in it within 365 days from the date of acquisition by the first acquirer.

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Mário Sérgio De Bortoli

Responsável Legal / Legal Representative / Responsable Legal / Responsabile Legale