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**Nome comercial:** FRESA




**EN Technical name:** Dental Surgical Drills / **Trade Name:** DRILL

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**EN** \* For all the types listed above, there is the conventional model and the irrigated model.

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**EN** Note: the difference between the conventional and irrigated drills is that the irrigated drills have a hole in them, as shown below:

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Nota: **ASTM F 899** - Standard Specification for Wrought Stainless Steels for Surgical Instruments

**EN** **Material used in the manufacture:** Stainless Steel ASTM F899.

Note: **ASTM F 899** - Standard Specification for Wrought Stainless Steels for Surgical Instruments

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
**EN** Enlarged image of the drills named “long”:

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**EN** The drills named “long” are exactly identical to the conventional drills, but the shanks are 5 mm longer.

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**EN** Lengths available for each diameter 5, 7, 9, 11, 13 and 15 mm.

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- 01 Fresa

### **EN 1- COMMERCIAL PRESENTATION FORMS OF THE PRODUCT**

The Implacil Drills are supplied in individual packs.

The primary packaging consists of a blister produced with virgin PET plastic, and it is sealed with Tyvek.

List of items contained inside the package

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### **EN 2- INDICATION OF USE AND PURPOSE**

The drills are used in dental surgeries in general.

The drills are tools that have a cutting edge, which through rotating movements (controlled by the counter-angle or manually), perforate the bone structure of the patient.

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**EN Helical / Irrigated Helical / Long Helical / Irrigated Long Helical:** Drill the bone and allow the placement of cylindrical implants.

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**EN Pilot / Irrigated Pilot / Long Pilot / Long Irrigated Pilot:** The pilot drills are used as intermediates between the helical drills, serving as a guide for subsequent drillings.

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**EN Tapered / Irrigated Tapered:** Drill the bone and allow the placement of tapered and morse taper implants.

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**EN Super Cutting / Irrigated Super Cutting / Long Super Cutting / Irrigated Long Super Cutting:** Serves to increase the diameter of the initial perforation in small dimensions.

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**EN Counter Sink / Irrigated Counter Sink / Long Counter Sink / Irrigated Long Counter Sink:** used to increase the diameter of the primitive drill, enabling the seating of the platform of the implant.

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**EN Reamer Drill / Irrigated Reamer Drill / Long Reamer / Irrigated Long Reamer:** Used to expand the bone bed for placement of the implants in a less traumatic way.

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**EN Spear / Irrigated Spear / Long Spear / Irrigated Long Spear:** serves to mark the site where the implant is to be inserted, promoting decorticalization in the bone crest.

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**EN Ball / Irrigated Ball / Long Ball / Irrigated Long Ball:** This drill is used for removing gingival or mucosal tissue.

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**EN Trephine / Irrigated Trephine / Long Trephine / Irrigated Long Trephine:** This drill is used for removing bone tissue by excision in fragments corresponding to its diameter.

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**EN Short Helical Guide Drill / Long Helical Guide Drill:** Drill the bone and allow the placement of cylindrical implants, used for guided surgeries.

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**EN Tapered Guide Drill:** Drill the bone and allow the placement of tapered and morse taper implants, used for guided surgeries.

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**EN Reamer Motor:** Used to expand the bone bed for placement of the implants in a less traumatic way.

The "irrigated" drills serve to prevent overheating of the bone in order to avoid its necrosis;

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**EN** The "long" drills are simply a length variation of the shank in relation to the conventional drills and their use remains at the discretion of the professional, because its choice will always depend on the space available for implant placement. For the surgical procedure, remove the drill from the instrument box, adapt it to the counter-angle through the shank, and check the lock of the drill on the counter-angle. With the drill in movement (200 to 600 RPM) start to drill. After reaching the desired depth, remove the counter-angle-Drill assembly from the bone cavity. Unlock the Drill and remove it from the counter-angle.

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**EN** They should be previously hygienized and sterilized according to the techniques below:

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**EN** The cleaning of a surgical instrument that has been in contact with tissues must be conducted in such a way that all parts of the item are exposed to the cleaning process.

To clean Implacil instruments, use only the following manual cleaning process.

Automated cleaning has not been validated. Do not use automated cleaning for Implacil instruments.

1. Initial rinse with cold water. The initial cold rinse is carried out in order to remove any gross contamination present in the item. Cold water is used, as temperatures above 60°C / 140 ° F will coagulate proteins, making it difficult to remove them from contaminated items.

2. Prepare the enzymatic detergent according to the recommendation \*.

3. Dip the products in the prepared detergent solution and let it soak for at least 5 minutes. While soaking, brush the items with a soft bristle brush.

4. Rinse: after removing the instruments from the detergent solution, rinse thoroughly for at least 3 times (for at least 1 minute) under running water.

5. Visually inspect each product for visible dirt. If any visible dirt remains, repeat the cleaning process (steps 3 and 4) until no visible dirt is left.

\* Recommendation for the cleaning step

Always use distilled, deionized or demineralized water for this procedure.

Never use saline solutions, especially sodium hypochlorite and saline, disinfectants, hydrogen peroxide or alcohol for cleaning or rinsing the cutters.

Never use steel wool or sponges and abrasive products, so that the cutters are not damaged.

Cleaning: Enzyme Detergent Composition: Protease, Lipase and Amylase enzymes

Dilution: 1mL for each liter of water;

Action time: 5 min (Manual);

Water temperature: 40 ° C to 55 ° C (Manual);

Note: The cleaning process must be conducted in such a way that all parts of the surgical instrument are exposed. Contamination remaining on an item after cleaning reduces the effectiveness of any subsequent decontamination or sterilization process.

Clean instruments as soon as possible after use. Do not allow blood and debris to dry on the instruments. If cleaning needs to be delayed, place groups of instruments in a container covered with enzyme solution or appropriate detergent to delay drying.

The instruments must be replaced as soon as the professional realizes their diminished effectiveness or degradation of the item.

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**EN** It is a procedure that aims at the total elimination of microorganisms, such as: viruses, bacteria, microbes, fungi (either in vegetative or sporulated form).

After the drying process, the instruments are sent to sterilization in an autoclave, using temperature and exposure time, according to the following method.

This product must be properly sanitized and sterilized before use. Please use for steam sterilization only the steam sterilization according to the following parameters:

Sterilization Parameters	
Temperature	132°C
Sterilization time	20 min
Drying time	20 - 60 min

Surgical instruments must be placed inside the container (perforated in the lid and at the bottom), or tray or box; provided that they can be prepared for the sterilization process. DO NOT open the autoclave to avoid rapid condensation; that is, let all the steam out first, making the Drying Cycle complete by itself. Thoroughly and periodically clean the autoclave, removing dirt and excess rust (iron oxide).

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**EN 5- EXPIRY DATE:** 5 years as from the date of manufacture.

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**EN 6- STORAGE CONDITIONS:** Keep at room temperature, in a cool place, away from sunlight and heat. Avoid exposure to heat and excess humidity. The drills should be stored clean, decontaminated and completely dry. The packaging in which the items are sterilized can provide an effective barrier to maintain the sterility of the item. After use, they should be kept in a clean, dry and cool place.

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**EN 7- TRANSPORT CONDITIONS:** The drills should be transported in a clean and dry place at room temperature. The transport should be performed appropriately to avoid falls and damages in their original pack.

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**EN 8- MANIPULATION CONDITIONS:** Maximum care must be taken in the cleaning, sterilization, accommodation (storage), transport and handling processes of the product to avoid mechanical shocks that change their anatomy and compromise their function.

Check the functioning of the product before initiating surgery. It should not be used if it has anatomical changes.

When handling the instruments (drills), use extreme caution in order to avoid injuries.

Unless indicated otherwise, sets of instruments are not sterile and should be sterilized before use.

Instruments should receive autoclave inside the instrument box according to 4.2.3b. Unpackaged instrument boxes (drills) do not maintain the sterility.

The drills should be replaced as soon as the professional notices its reduced efficacy.

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**EN 9- WARNINGS:** Exclusive use by the dental professional.

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**EN 10- PRECAUTIONS:** Incorrect maintenance, cleaning and maintenance can make the instruments inappropriate for their intended use, cause corrosion, dismounting, distortion and/or break or even cause injuries on the patient or operating team.

Below is a list of potential complications:

- Injury or pain in the soft tissues, in case of incorrect use or breakage of the instruments.
- Infection, if the instrument is not conveniently clean and sterilized.
- Damages to the nerves or nearby organs as a result of slide or bad positioning of a damaged instrument.
- Damages caused by instruments used to screw, due to excess forces that occur during use.
- Increase in the time of operation due to dismounting of the instruments during surgery.
- Due to the necessary mechanical functionalities, majority of the instruments (drills) are made in non-implantable materials. In case an instrument breaks, no fragment should remain inside the patient because it can cause post-operative complications, such as allergies, infections or complications of the biological kind, associated to the release of metallic components, possibly requiring future intervention.

To avoid environmental contamination or undue use of the Implacil non-articulated cutting surgical instruments, when these are unusable, the instruments should be disposed in the appropriate place, i.e., in a place for potentially contaminated materials.

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