**PORTUGUÊS**

Pilar Protético Ideale

**INGLÊS**

Ideale Abutment Prosthetic

**ESPAÑHOL**

Pilar Protético Ideale

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Revisão: 04

## PRODUCT:

Technical name: Prosthetic Abutment

Commercial name: Ideale Abutment Prosthetic

## ATTENTION!!!

Single-use product. It is recommended not to reuse.

Do not use if packaging is broken or damaged.

Non-Sterile Product. Sterilize before use.

Reprocessing prohibited.

## Commercial models:

IDEALE ABUTMENT CM

IDEALE ABUTMENT CM AR

IDEALE ABUTMENT SCREW Through

IDEALE ABUTMENT ANGLE CM

IDEALE ABUTMENT ANGLE CM AR

COPING IDEALE MORSE TAPER

IDEALE ABUTMENT HI

IDEALE ABUTMENT ANGLE HI

IDEALE ABUTMENT HE

IDEALE ABUTMENT ANGLE HE

## 4. INDICATION OF USE AND PURPOSE

The abutments are indicated for the final manufacture of a single or multiple prosthesis element, where there is a need for functional and/or aesthetic rehabilitation, being fixed using a screw or specific cement for dental procedures and being able to make the necessary corrections in the relationship of the implants to the axis of the teeth, both in the vertical dimension and in angulation.

## 5. HOW TO USE: DETAILED DESCRIPTION

### 2.1 IDEALE ABUTMENT CM

Consisting of a abutment straight conical, manufactured in commercially pure titanium grade 4, according to ASTM-F67 standard.

The abutments are for single use, presented with strap heights of 0.0mm; 0.8mm; 1.5mm; 2.5mm; 3.5mm; 4.5mm; 5.5mm.

Its upper base has a thread to screw the prosthesis in. The lower base has an M1.8 thread to adapt to the implant.


It is used in the final production of a screwed prosthesis element, whether single or multiple.

The product is individually packaged, each package consists of an Ideale abutment cone morse (available in gold color).

Installation torque: 30 N.cm

Compatibility : Implant UNIII Cone Morse AR

### 2.2 IDEALE ABUTMENT CM AR

Consisting of a abutment straight conical, manufactured in commercially pure titanium grade 4, according to ASTM-F67 standard.

And a screw for fixing the abutment to the implant, made of ASTM-F136 titanium alloy.

The abutments are for single use, presented with strap heights of 0.0mm; 0.8mm; 1.5mm; 2.5mm; 3.5mm; 4.5mm; 5.5mm.

At its upper base it has a platform for fitting the prosthesis and at its lower base it has a hexagon with a positioning function and has a through screw with M1.8 thread for fixing the abutment to the implant.

It is used in the final production of a screwed prosthesis element, whether single or multiple.

The product is individually packaged, each package consists of an ideale abutment cone morse AR (available in gold color) and a fixing screw.

Installation torque: 20 N.cm

Compatibility : Implant UNIII Cone Morse AR

### 2.3 IDEALE ABUTMENT SCREW THROUGH

Consisting of a abutment straight conical, manufactured in commercially pure titanium grade 4, according to ASTM-F67 standard.

And a screw for fixing the abutment to the implant, made of ASTM-F136 titanium alloy.

The abutments are for single use and are presented with strap heights of 0.8mm; 1.5mm; 2.5mm; 3.5mm; 4.5mm and 5.5mm.

This abutment has a through screw (M1.8 thread), on its upper base it has an M1.4 thread for fitting the prosthesis and on the lower base it has a screw through with M1.8 thread for fixing the abutment to the implant.

It is used in the final production of a screwed prosthesis element, whether single or multiple.

The product is individually packaged, each package consists of an Ideale Abutment Screw Through (available in gold color) and a fixing screw.

Installation torque: 20 N.cm

Compatibility : Implant UNIII Cone Morse AR

## 2.4 IDEALE ABUTMENT ANGLE CM

Consisting of an abutment angled conical, made of ASTM-F136 titanium alloy.

And a screw for fixing the abutment to the implant, made of ASTM-F136 titanium alloy.

The abutments are for single use and are presented with belt heights of 1.5mm, 2.5mm and 3.5mm. They have an angle of 17° or 30°.

Its upper base has an M1.4 thread for screwing the prosthesis. Its lower base has a through screw with an M1.8 thread for fixing the abutment to the implant.

It is used in the final production of a screwed prosthesis element, whether single or multiple.

The product is individually packaged, each package consists of an Ideale abutment Angle CM (available in gold) and a fixing screw.

Installation torque: 20 N.cm

Compatibility : Implant UNIII Cone Morse AR

## 2.5 IDEALE ABUTMENT ANGLE CM AR

Consisting of an angled conical abutment, made of ASTM-F136 titanium alloy, and a fixing screw made of ASTM-F136 titanium alloy.

The abutments are for single use and are presented with belt heights of 1.5mm, 2.5mm and 3.5mm. They have an angle of 17° or 30°.

Its upper base has an M1.4 thread for screwing the prosthesis. Its lower base has a hexagon with a positioning function and has a through screw with an M1.8 thread for fixing the abutment to the implant.

It is used in the final production of a screwed prosthesis element, whether single or multiple.

The product is individually packaged, each package consists of an Ideale Abutment Angle CM AR (available in gold color) and a fixing screw.

Installation torque: 20 N.cm

Compatibility : Implant UNIII Cone Morse AR

## 2.6 COPING IDEALE MORSE TAPER

Consisting of a straight cylindrical abutment, made of commercially pure titanium grade 4, according to ASTM-F67 standard, and a fixing screw made of titanium alloy ASTM-F136.

The coping are for single use only. The upper base has a hole through which the screw passes in order to fix the hood to the abutment. The lower base has an internal fitting, which is used to adapt to the Ideale Abutment of the Implacil System.

Its height can be 4 or 6mm and the base diameter can be 3.3 or 4.5


It is used in the final production of a screwed prosthesis element, whether single or multiple.

The product is individually packaged, each package consists of a titanium abutment coping ideale morse cone AR (available in gold color) and a titanium alloy fixing screw.

Installation torque: 20 N.cm

Compatibility: Ideale Abutment Cone Morse; Ideale Abutment Cone Morse AR; Ideale Abutment Screw Through; Ideale Abutment Angle Cone Morse; Ideale Abutment Angle Cone Morse AR

## 2.7 IDEALE ABUTMENT HI

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Consisting of a straight conical abutment, made of ASTM-F136 titanium alloy and a fixing screw made of ASTM-F136 titanium alloy.

The abutments are for single use, presented with strap heights of 1.0mm; 2.0mm; 3.0mm; 4.0mm; 5.0mm.


It is used in the final production of a screwed prosthesis element, whether single or multiple.

The product is individually packaged, each package consists of an ideale abutment HI (available in gold color) and a fixing screw.

Installation torque: 30 N.cm

Compatibility: Implant connection HI

## 2.8 IDEALE ABUTMENT ANGLE HI

Consisting of an abutment angle conical, made of ASTM-F136 titanium alloy and a fixing screw made of ASTM-F136 titanium alloy.

The pillars are for single use and are presented with strap heights of 1.0mm, 2.0mm and 3.0mm.

They have an angle of 17° or 30°.

Its upper base has an M1.4 thread for screwing the prosthesis. Its lower base has a through screw with an M1.8 thread for fixing the abutment to the implant.

It is used in the final production of a screwed prosthesis element, whether single or multiple.

The product is individually packaged, each package consists of an ideale abutment angle HI (available in gold color) and a fixing screw.

Installation torque: 20 N.cm

Compatibility: Implant connection HI

## 2.9 IDEALE ABUTMENT HE

Consisting of a abutment straight conical, made of ASTM-F136 titanium alloy and a fixing screw made of ASTM-F136 titanium alloy.

The abutments are for single use, presented with strap heights of 1.0mm; 2.0mm; 3.0mm; 4.0mm; 5.0mm.


It is used in the final production of a screwed prosthesis element, whether single or multiple.

The product is individually packaged, each package consists of an ideale abutment HE (available in gold) and a fixing screw.

Installation torque: 30 N.cm

Compatibility: Implant connection HE

## 2.10 IDEALE ABUTMENT ANGLE HE

Consisting of an abutment angle conical, made of ASTM-F136 titanium alloy and a fixing screw made of ASTM-F136 titanium alloy.

The abutments are for single use and are presented with strap heights of 1.0mm, 2.0mm and 3.0mm.

They have an angle of 17° or 30°.

Its upper base has an M1.4 thread for screwing the prosthesis. Its lower base has a through screw with an M1.8 thread for fixing the abutment to the implant.

It is used in the final production of a screwed prosthesis element, whether single or multiple.

The product is individually packaged, each package consists of an ideale abutment Angle HE (available in gold color) and a fixing screw.

Installation torque: 20 N.cm

Compatibility: Implant connection HE

Note: See Annex I to consult COMMERCIAL MODELS AND REFERENCE CODES

## 3.FORMS OF COMMERCIAL PRESENTATION OF THE PRODUCT

The Ideale Abutment Prosthetic Implacil De Bortoli is presented in individual, non-sterile packaging.

The primary packaging consists of a sealed aluminum sachet, and its final packaging is a cardboard envelope, which allows the internal contents to be viewed.

The product is presented in individual packaging. Package contents:

01 Ideale Abutment Prosthetic

01 Fixing Screw (when applicable)

03 Traceability Labels

All models come with 3 traceability labels, which must be attached: to the clinical record, to the document to be given to the patient, and to the tax documentation that generates the charge.

Each traceability label contains:

- a) commercial name or model;
- b) identification of the manufacturer;
- c) product or system component code;
- d) batch number and ANVISA registration number.

This instruction for use will be available on the Implacil De Bortoli website ([www.implacil.com.br](http://www.implacil.com.br)) in accordance with current legislation on the provision of instructions for use in non-printed format for health products – RDC No. 751 of September 15, 2022. The printed instruction for use can be requested through the Customer Service – 0800.779.0011, and will be sent at no additional cost.

To locate it, check the “IFU” number and the version of the instruction that is indicated on the Label of the purchased product. To read the instructions for use on the website, you need the Adobe Acrobat Reader program.

User Alert: Always observe the correlation of the version of the instructions for use indicated on the label of the purchased product, as made available by Implacil De Bortoli.

#### Packaging Type:

The material and type of packaging used to package the Ideale Abutments Prosthetic Implacil De Bortoli are well known and universally used, both for medicines and health products, and are listed in the ANVISA Packaging Table.

#### Packaging features and dimensions:

Primary Packaging: Blister sealed with Tyvek 50mm X 80mm  
ANVISA Packaging Table: PVC Blister – code 0119 Tyvek ®  
Envelope – code 0332

Secondary packaging: Carton packaging - 55 mm x 85 mm

#### **4. PHYSICAL PRINCIPLE AND FUNDAMENTALS OF THE PRODUCT'S TECHNOLOGY, APPLIED TO ITS OPERATION AND ACTION:**

Abutments Prosthetic are mechanical components that act as a link between the implant and the dental prosthesis, providing support for the latter and making corrections to the height and angle between the implant and the tooth. They are fixed to the implant using a thread. The prosthetic elements can be fixed to the abutments using cement or screws, and it is up to the professional to decide which method is best to use.

#### **5. OPTIONALS, CONSUMABLES AND SUPPORT MATERIALS**

Optional: Not applicable

Consumables: Not applicable

Support materials: Instructions for Use and traceability labels.

Manufacturing items for exclusive use in the company's products:

All parts described in these Instructions for Use are for exclusive use in the Prosthetic Abutment product.

Attention: The use of any part, accessory or material not specified or provided for in these Instructions for Use is the sole responsibility of the user.

#### **6. INSTRUCTIONS FOR USE OF THE PRODUCT**

- Abutments Prosthetic Implacil come packaged in a blister sealed with Tyvek®.
- Check that the packaging is in perfect condition. Reject the product if the packaging is damaged.
- Open the blister and remove the abutment from inside the blister, check that it is perfectly sealed.
- The product is not sterile, and this procedure must be carried out by a professional in accordance with current standards (NBR 14332), or any current legislation or standard that may complement or replace it.
- Screw the abutment into the implant.
- Check the model description (item C of these Instructions for Use) for the maximum torque to be applied.
- The Abutment has been previously sanitized in a classified room. If you wish to sanitize it again, follow these steps:
  - Initial rinse of the item with cold water.
  - Wash with a neutral to slightly basic detergent
  - Final rinse with deionized water or water purified by reverse osmosis.
- The product is not sterile, this procedure must be carried out by a professional.
- Sterilization is a procedure that aims to completely eliminate microorganisms, such as viruses, bacteria, microbes, and fungi (whether in vegetative or sporulated form). After the drying process, the abutments must be sterilized in an autoclave, using the temperature and exposure time according to the following method: Sterilization temperature: 132°C, Sterilization cycle: 20 minutes, and Drying cycle: 20 – 60 minutes. The abutments must be placed inside the container (perforated in the lid and bottom), or on a tray or in the box; provided that they are arranged in such a way that they can be prepared for the sterilization process. DO NOT open the autoclave to avoid rapid condensation; that is, let all the steam escape first, allowing the Drying Cycle to complete itself. Clean the autoclave thoroughly and periodically, removing dirt and excess rust (iron oxide).

It is the dentist's responsibility to sterilize the product as indicated in these instructions for use. Failure to sterilize the product properly may cause biocontamination with infectious diseases, resulting in failure to use the product.

##### **6.1 TRACEABILITY LABEL:**

Abutments Prosthetic Implacil come with three traceability labels that must be attached 01 to the medical record, 01 to

the tax document to be given to the patient and 01 to the tax collection document.

Traceability is done through the product batch.

## **7. COMPOSITION**

As raw material for the Abutments, the following are used: Titanium Alloy according to ASTM F136 standard

## **8. VALIDITY PERIOD**

Indefinite validity

## **9. SPECIAL CONDITIONS FOR TRANSPORTATION, STORAGE AND PACKAGING OF THE PRODUCT:**

Store at a temperature between 15°C and 45°C and humidity between 35% and 85%, protected from light and heat. Avoid exposure to heat and excessive humidity. Do not place excessive load on the packaging to avoid compromising the packaging seals.

Abutments Prosthetic Implacil De Bortoli must be transported in a clean, dry place at room temperature. Transportation must be carried out appropriately to avoid falling and damaging the original packaging and mechanical shocks that may damage the product.

### **ATTENTION!!!**

Single-use product. It is recommended not to reuse.

Do not use if packaging is broken or damaged.

Non-Sterile Product. Sterilize before use.

Reprocessing prohibited.

## **10. WARNINGS AND/OR PRECAUTIONS TO BE ADOPTED**

Abutments Prosthetic Implacil should only be used by professionals specialized in the area. It is necessary to develop a prosthetic execution plan, observing the diameter, angle and length for the adaptation of the prosthetic component to the implant in order to avoid future complications.

Abutments Prosthetic Implacil are designed for single use and cannot be resterilized.

Reuse of this product may cause adverse biological effects resulting from product residues, microorganisms and/or substances resulting from previous use and/or reprocessing, changes in the original physical, mechanical and chemical, macro and micro structural characteristics of the product that may impair its intended functionality. Reuse of this product does not guarantee its safety and effectiveness and exempts any warranty of related products.

Do not use it under any circumstances if the packaging is broken. The product is not sterile, and this procedure must be carried out by a professional before any intervention on the patient.

The professional must observe the force exerted when applying the product so as not to damage it and inform the patient of the need for periodic monitoring to avoid physical

and mechanical efforts in order to not subject the product to undue stress, as well as the need for correct cleaning.

## **11. UNDESIRABLE SIDE EFFECTS AND CONTRAINDICATIONS**

The final result of the application of the product may not be successful in cases where the technique used is not appropriate and the patient does not undergo the indicated examinations and proper prosthetic planning, which may result in loss or fracture of the product. Incorrect application of the product may cause side effects in the region where it was applied, such as: pain, swelling, sensitivity, tissue reaction and infection.

It is the professional's responsibility to advise the patient on the need for professional monitoring after the procedure and to follow the guidelines on care, hygiene and medication prescription.

The use of Abutments Prosthetic Implacil is contraindicated in cases where the patient presents chronic periodontal inflammation, inadequate oral hygiene and parafunctional habits, active intraoral infection and in cases of Immediate Loading in which the primary stability of the implant is not adequate.

## **12. PRODUCT SAFETY AND EFFICACY:**

Abutments Prosthetic Implacil are completely safe as long as safety rules are followed in all procedures performed. The product was developed with non-toxic medical -dental material. When stored and preserved according to the instructions for use, the product will not undergo changes or modifications in its characteristics.

## **13. INSTALLATION OR CONNECTION TO OTHER PRODUCTS**

Abutments Prosthetic Implacil are connected to previously installed and osseointegrated implants, according to the compatible models for each type of connection.

The use of any part, accessory or material not specified or provided for in these instructions for use is the sole responsibility of the user.

## **14. PRECAUTIONS IN CASE OF DISPOSAL**

To avoid environmental contamination or improper use of Abutments Prosthetic, when they are no longer usable, they must be disposed of in hospital waste or in a location for potentially contaminated materials.

## **15. WARRANTY TERMS**

IMPLACIL DE BORTOLI - MATERIAL ODONTOLÓGICO SA guarantees the owner of this product against any material or manufacturing defect that may arise within 365 days from the date of purchase by the first purchaser. If transported, stored and packaged properly in accordance with these Instructions for Use, the product will have an indefinite shelf life.

It is estimated that dental implants have a useful life of approximately 15 years, therefore, as determined by Law No. 8078 of September 11, 1990, as long as there is a complaint within the ninety-day statute of limitations after the defect is evidenced, even if the defect only manifests itself after the end of the contractual warranty period , Implacil De Bortoli undertakes to replace the device with an identical or similar device.

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© Diámetro Plataforma