**PORTUGUÊS**

FRESA

**INGLÊS**

DRILL

**ESPAÑHOL**

FRESA

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Revisão: 09

## PRODUCT

Technical Name: Dental Surgical Drills

Commercial Name: DRILL

## ATTENTION!!!

**Do not use if the packaging is broken or damaged.**

**Non-Sterile Product. Sterilization is required before use.**

**Product subject to reprocessing.**

## Commercial Models:

MALE INTERNAL THREAD

MALE INTERNAL THREAD M2.0/M1.8/M1.6

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MALE EXTERNAL THREAD

DRILL 1.2

34583 - MALE EXTERNAL THREAD STRETTO IMPLANT

27212 - MINI DRILL SPEAR 2.0X05

25140 - MINI DRILL SPEAR 2.0X06

25157 - MINI DRILL SPEAR 2.0X07

25164 - MINI DRILL SPEAR 2.0X08

25027 - MINI DRILL SPEAR 2.0X10

## 1- INDICATION OF USE AND PURPOSE

Drills are used in orthodontic surgeries in general.

## 2- HOW TO USE: DETAILED DESCRIPTION

MALE INTERNAL THREAD

Diameter: Ø2mm. Length: 22mm

MALE INTERNAL THREADED DIGITAL

Diameter: Ø2mm. Length: 28mm

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Diameter: Ø 1.4; Ø2.0mm. Lengths 5.0; 6.0; 7.0; 8.0; 10.0mm.

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Diameters: Ø 8.0 and Ø 10.0mm. Lengths: 6.0; 7.0; 8.0; 10.0mm.

MALE EXTERNAL THREAD

Diameter: Ø3.3; Ø3.75; Ø4.0; Ø4.75mm . Length: 26mm

DRILL 1.2

Diameter: Ø1.2mm. Length: 26mm

34583 - MALE EXTERNAL THREAD STRETTO IMPLANT

Diameter: Ø3.2mm. Length: 26.5mm

For the surgical procedure, remove the bur from the instrument box, adapt it to the contra-angle using the rod, and check that the drill bit is locked in the contra-angle. With the drill bit in motion (200 to 600 RPM), start drilling. After reaching the desired depth, remove the contra-angle-burr assembly from the bone cavity. Unlock the burr and remove it from the contra-angle .

They must be previously cleaned and sterilized using the following techniques:

### a) Cleaning:

Cleaning of a surgical instrument (burr) that has been in contact with tissue must be carried out in such a way that all parts of the item are exposed to the cleaning process.

The cleaning process should involve an initial rinse of the item with cold water, a wash with a mild to slightly basic, low-sudsing detergent, a rinse with water (tap or conditioned) followed by a final rinse with deionized or reverse osmosis purified water.

The initial cold rinse is conducted in order to remove any gross contamination present on the item.

Cold water is used because temperatures above 140°F will coagulate proteins, making them difficult to remove from contaminated items.

After rinsing in cold water, items should be cleaned with a neutral to slightly basic detergent with low foaming and a pH between 7.0 and 10.0.

If the contamination contains heavy organic residues (protein or fat) or dried dirt, an enzymatic detergent can be used.

The detergent used should also be easy to rinse and not leave any residual compounds. A low-foaming detergent reduces the formation of aerosol droplets, which can carry microorganisms. Detergent residues can cause stains and interfere with the action of chemical disinfectants.

Washed items should be rinsed to remove detergent, preferably with a final rinse in deionized or reverse osmosis purified water.

The cleaning process should be conducted in such a way that all parts of the surgical instrument (bur) are exposed. Contamination remaining on an item after cleaning reduces

the effectiveness of any subsequent decontamination or sterilization process.

Clean instruments (burs) as soon as possible after use. Do not allow blood and debris to dry on instruments (burs). If cleaning must be delayed, place groups of instruments (burs) in a covered container with enzymatic solution or appropriate detergent to delay drying.

#### b) Sterilization

It is a procedure that aims to completely eliminate microorganisms, such as viruses, bacteria, microbes, fungi (whether in vegetative or sporulated form). After the drying process, the instruments (burs) are sent for sterilization in an autoclave, using temperature and exposure time, according to the following method:

Sterilization temperature: 132°C, Sterilization cycle: 20 minutes and Drying cycle: 20 – 60 minutes.

Surgical instruments must be placed inside the container (perforated in the lid and bottom), or on a tray or in the box; provided they are arranged in such a way that they can be prepared for the sterilization process.

DO NOT open the autoclave to avoid rapid condensation; that is, let all the steam escape first, allowing the Drying Cycle to complete itself.

Clean the autoclave thoroughly and periodically, removing dirt and excess rust (iron oxide)

### **3. FORMS OF COMMERCIAL PRESENTATION OF THE PRODUCT**

Implacil Strawberries are presented in individual packaging. The primary packaging consists of a blister made from virgin PET plastic and is sealed with Tyvek.

List of items contained inside the package:

- 01 Drill

This instruction for use will be available on the Implacil De Bortoli website ([www.implacil.com.br](http://www.implacil.com.br)) in accordance with current legislation on the provision of instructions for use in non-printed format for health products – RDC 751/2022. The printed instruction for use will be available to the customer who requests it, and will be sent free of charge. Printed instructions for use can be requested through Customer Service – 0800.779.0011, and will be sent at no additional cost.

To locate it, check the “IFU” number and the version of the instruction that is indicated on the Label of the purchased product. To read the instructions for use on the website, you need the Adobe Acrobat Reader program.

User Alert: Always observe the correlation of the version of the instructions for use indicated on the label of the purchased product, as made available by Implacil De Bortoli.

### **4. PHYSICAL PRINCIPLE AND FUNDAMENTALS OF THE PRODUCT'S TECHNOLOGY, APPLIED TO ITS OPERATION AND ACTION:**

Burs are tools that have a cutting edge, which through a rotary movement (controlled by the contra-angle or manually), causes the drilling of the patient's bone structure.

### **5. COMPOSITION**

As raw material for the Drills, ASTM 899 Stainless Steel is used.

### **6. VALIDITY PERIOD**

Undetermined

### **7. SPECIAL CONDITIONS FOR TRANSPORTATION, STORAGE AND HANDLING OF THE PRODUCT:**

Store at room temperature, in a cool place, away from light and heat. Avoid exposure to heat and excessive humidity.

Straws should be stored clean, decontaminated and completely dry. The packaging in which the items are sterilized can provide an effective barrier to maintain the sterility of the item. After use, they should be stored in a clean, dry and cool place.

Strawberries must be transported in a clean, dry place at room temperature. They must be transported properly to avoid falling and damaging their original packaging.

The utmost care must be taken in the process of cleaning, sterilization, storage, transportation and handling of the product to avoid mechanical shocks that alter its anatomy and impair its functionality.

Before starting surgery, check that the product is working properly. It should not be used if there are any anatomical changes.

When handling instruments (burs), use extreme caution to avoid injury.

Unless otherwise indicated, instrument sets are non-sterile and must be sterilized prior to use.

Instruments must be autoclaved inside the instrument box. Instrument boxes (burs) when not packaged do not maintain sterility.

The cutters must be replaced as soon as the professional notices their reduced effectiveness.

Exclusive use by dental professionals.

### **8. WARNINGS AND/OR PRECAUTIONS TO BE ADOPTED**

Exclusive use by dental professionals.

Incorrect maintenance, cleaning and servicing may render the instruments unsuitable for their intended use, cause corrosion, disassembly, distortion and/or breakage or even cause injury to the patient or operating team.

Below is a list of potential complications:

- Injury or pain in soft tissues, in case of incorrect use or breakage of instruments.

- Infection, if the instruments are not properly cleaned and sterilized.
- Damage to nearby nerves or organs as a result of slipping or mispositioning of a faulty instrument.
- Damage caused by the instruments used for screwing, due to excessive forces that occur during use.
- Increased operating time due to having to disassemble the instruments during surgery.
- Due to the necessary mechanical functions, most instruments (burs) are made of non-implantable materials. In the event of an instrument breaking, no fragments should remain in the patient, as this may cause post-operative complications, such as allergies, infections, or complications of a biological nature, associated with the release of metal components, possibly requiring subsequent intervention. To avoid environmental contamination or improper use of Implacil non-articulated cutting instruments, when they are to be discarded, they must be discarded in an appropriate place, that is, in a place for potentially contaminated materials.

#### **9. CONTRAINDICATIONS**

Not applicable

#### **10. ADVERSE EFFECTS**

Not applicable

#### **11. PRECAUTIONS IN CASE OF DISPOSAL**

To avoid environmental contamination or misuse of Prosthetic Abutments, when they are no longer used, they must be disposed of in hospital waste or in a location for potentially contaminated materials.

#### **12. WARRANTY TERMS**

IMPLACIL DE BORTOLI - MATERIAL ODONTOLOGICO SA guarantees the owner of this product a guarantee against any material or manufacturing defect that may arise within 365 days, counting from the date of purchase by the first purchaser.

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